

Certificate in relation to air transportation by passengers with walking disabilities

We hereby certify that our product

Twion M24 / T24

including the built-in lithium-ion battery fulfil the requirements of

IATA chapter 2.3.2.4 (Wheelchair / Mobility Aids)

The lithium-ion batteries are tested according to UN Manual of Tests and Criteria, Part III, Subsection 38.3.

The system may be carried as checked baggage in the cargo compartments.

For lithium batteries installed in mobility aids, there is no Watt-hour limit.

Additional information: The rated energy of the built-in battery (not removable by the user) is < **160 Wh**.

Condition precedent to meet the above listed requirement is the transportation of the twion wheels in a way that an unintentionally switching on and damage by other luggage is avoided. We recommend using the dedicated Alber transport axle.

We strongly recommend contacting the airline sufficiently in advance to request approval for transport.

This certificate does not apply to damaged or defective batteries.

Due to the frequent amendments of the respective regulations, which partly occur on an annual basis this certificate is valid until **31st December 2025**.

Albstadt, 22nd of November 2024

Alber GmbH



Johannes Beltman
Managing Director